

# Deming Designs Beach Walkers

1090 Cobblestone Dr  
Pensacola, FL 32514

850-478-5765  
www.beachwheelchair.com

## Stainless Steel Beach Walker



### Specifications

Frame:	316L Stainless Steel Tubing
Wheels	30cm Wheeliez Wheels
Casters	Rear are fixed, Front rotate 360 degrees
Weight	15lbs to 35lbs
Width	30" outside to outside
Height	Standard 32" can be adjustable. Children sizes are also available
Weight Capacity	Standard 350lbs


## Aluminum Beach Walker



### Specifications

Frame	Powder Coated Aluminum
Wheels	30cm Wheeliez Wheels, in Front 22cm Wheeliez Wheels and 20cm Polycarbonate Wheels in Rear
Casters	316l Stainless Steel
Weight	28lbs
Width	Width clearance 33 inches
Height	36.5 inches to 41 inches
Weight Capacity	250lbs

## Standing Chart and Pricing



Product		Price	Subtotal
<b>Aluminum Walker</b> <small>Please indicate color choice no additional charge</small>	<b>Red Blue</b>	<b>\$895.00</b>	
<b>Stainless Steel Walker</b>		<b>\$995.00</b>	
<b>Quick Release Pin Option</b> <small>(Stainless Only)</small>		<b>\$90.00</b>	
<b>Custom Powder Coating</b> <small>(Stainless Only)</small>	<b>Your choice of color</b>	<b>\$300.00</b>	
<b>Umbrella with Stainless Steel Holder</b>		<b>\$85.00</b>	
<b>Frieght</b>			
<b>Total</b>			

## Billing and Shipping Information

### Billing

First Middle Last

Street State Zip

Phone Fax Email

### Shipping

First Middle Last

Street State Zip

Phone Fax Email