

DEBUG TILT-IN-SPACE RECLINING ALL-TERRAIN WHEELCHAIR



**Images shown features optional add-on seating accessories*

FOR MORE INFORMATION WWW.BEACHWHEELCHAIR.COM

DeBug Mobility Products
78 E Nine Mile Road
Pensacola, FL, 32534

850-478-5765
office@debugmobilityproducts.com
www.beachwheelchair.com

DEBUG TILT N RECLINE ALL-TERRAIN WHEELCHAIR

The DeBug Tilt N Recline All Terrain Wheelchair is engineered for individuals who require advanced postural support without sacrificing outdoor freedom. Built from high grade 316L stainless steel and designed with our unique rear axle suspension, this chair delivers unmatched stability and comfort across sand, soft terrain, and shoreline environments. With fully adjustable tilt in space and recline functions, caregivers can fine tune positioning for pressure relief, respiratory comfort, and long duration seating needs. Every detail reflects DeBug's commitment to safety, durability, and inclusive outdoor access. Empowering users to enjoy beaches, parks, and waterfront destinations with confidence and ease

- LONG LASTING #316L STAINLESS STEEL FRAME
- REMOVABLE TELESCOPIC FOOTREST
- FITS THROUGH A STANDARD ADA DOORWAY
- POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS GRASS, GRAVEL, SAND OR SNOW
- DUAL STITCHED SLING WITH TENERA 130 THREAD
- PHIFERTEX PLUS SLING MATERIAL IS MACHINE WASHABLE AND RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE

SPECIFICATIONS

Seat Sizes	ENGLISH (IN)	METRIC (CM)
Widths	16, 18, 22	41, 46, 56
Depth	17	43
Seat-to-Floor height	17, 21	43, 53
Backrest Height	18	46
Telescopic Footrest <i>Knee-to-Heel</i>	15-22	38-56
Weight Capacity	Up to 300 lbs	Up to 136 kg
Chair Width	33.5 inches	85 cm
Chair Weight	60-65 lbs	27-29 kg