

DEBUG TILT-IN-SPACE ALL-TERRAIN WHEELCHAIR



**Image shown features optional add-on seating accessories*

FOR MORE INFORMATION WWW.BEACHWHEELCHAIR.COM

DeBug Mobility Products
78 E Nine Mile Road
Pensacola, FL, 32534

850-478-5765
office@debugmobilityproducts.com
www.beachwheelchair.com

DEBUG TILT N SPACE

ALL-TERRAIN WHEELCHAIR

The DeBug Tilt-in-Space is built for comfort, stability, and advanced clinical support. It features a corrosion-resistant 316L stainless steel frame, and soft-terrain Wheelleez® balloon wheels. With an incremental tilt-in-space system for pressure relief, blood-pressure regulation, and improved torso support, it allows users to enjoy the beach and outdoors safely, comfortably, and with confidence. Whether it's a day at the beach, a lakeside outing, or navigating soft ground at parks and trails, the DeBug Beach Wheelchair delivers stability, comfort, and independence.

- LONG LASTING #316L STAINLESS STEEL FRAME
- REMOVABLE TELESCOPIC FOOTREST
- FITS THROUGH A STANDARD ADA DOORWAY
- POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS GRASS, GRAVEL, SAND AND SNOW
- DUAL STITCHED SLING WITH TENERA 130 THREAD
- PHIFERTEX PLUS SLING MATERIAL IS MACHINE WASHABLE AND RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE

SPECIFICATIONS

Seat Sizes	ENGLISH (IN)	METRIC (CM)
Widths	16, 18, 22	41, 46, 56
Depth	17	43
Seat-to-Floor height	17, 21	43, 53
Backrest Height	18	46
Telescopic Footrest <i>Knee-to-Heel</i>	15-22	38-56
Weight Capacity	Up to 300 lbs	Up to 136 kg
Chair Width	33.5 inches	85 cm
Chair Weight	60-65 lbs	27-29 kg