

DEBUG FIXED LEG REST AQUATIC WHEELCHAIR



FOR MORE INFORMATION WWW.BEACHWHEELCHAIR.COM

DeBug Mobility Products
78 E Nine Mile Road
Pensacola, FL, 32534

850-478-5765
office@debugmobilityproducts.com
www.beachwheelchair.com

DEBUG FIXED LEG REST AQUATIC WHEELCHAIR

The DeBug Aquatic Wheelchair is designed for use in pools, splash pads, waterparks, and shower facilities while also functioning as a comfortable manual wheelchair on paved surfaces. Constructed from 316L stainless steel, it helps facilities provide accessible water entry while meeting ADA accessibility guidelines. With its optional interchangeable beach conversion kit, the chair can be quickly transformed from an aquatic wheelchair into a beach wheelchair, giving users greater freedom in and around the water.

- LONG LASTING #316L STAINLESS STEEL FRAME
- REMOVABLE TELESCOPIC FOOTREST
- FITS THROUGH A STANDARD ADA DOORWAY
- LOW MAINTENANCE DELRIN BUSHINGS
- TUBELESS SOLID CORE TIRES, POLYCARBONATE MAG WHEELS
- DUAL STITCHED SLING WITH TENERA 130 THREAD
- PHIFERTEX PLUS SLING MATERIAL IS MACHINE WASHABLE AND RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE

SPECIFICATIONS

Seat Sizes	ENGLISH (IN)	METRIC (CM)
Widths	16, 18, 22	41, 46, 56
Depth	17	43
Seat-to-Floor height	17, 21	43, 53
Backrest Height	18	46
Telescopic Footrest <i>Knee-to-Heel</i>	15-22	38-56
Weight Capacity	Up to 300 lbs	Up to 136 kg
Chair Width	33.5 inches	85 cm
Chair Weight	55-60 lbs	25-27 kg