

DEBUG RECLINING ALL-TERRAIN WHEELCHAIR



FOR MORE INFORMATION WWW.BEACHWHEELCHAIR.COM

**Image shown may include optional accessories*

DeBug Mobility Products
78 E Nine Mile Road
Pensacola, FL, 32534

850-478-5765
office@debugmobilityproducts.com
www.beachwheelchair.com

DEBUG RECLINING ALL-TERRAIN WHEELCHAIR

Premium beach mobility built for comfort and durability. The 316L marine-grade stainless steel frame resists rust and salt, while Wheeleez® balloon wheels glide over soft sand with ease. The reclining backrest provides customizable support, improved circulation, and all-day comfort at the beach. Strong, stable, and designed for accessible shoreline enjoyment.

- LONG LASTING #316L STAINLESS STEEL FRAME
- REMOVABLE TELESCOPIC LEG REST
- FITS THROUGH A STANDARD ADA DOORWAY
- POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS GRASS, GRAVEL, SAND, AND SNOW.
- DUAL STITCHED SLING WITH TENERA 130 THREAD
- PHIFERTEX PLUS SLING MATERIAL IS MACHINE WASHABLE AND RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE

SPECIFICATIONS

Seat Sizes	ENGLISH (IN)	METRIC (CM)
Widths	16, 18, 22	41, 46, 56
Depth	17	43
Seat-to-Floor height	17, 21	43, 53
Backrest Height	18	46
Telescopic Footrest <i>Knee-to-Heel</i>	15-22	38-56
Weight Capacity	Up to 300 lbs	Up to 136 kg
Chair Width	33.5 inches	85 cm
Chair Weight	60-65lbs	27-29 kg