

DEBUG XL FRAME ALL-TERRAIN WHEELCHAIR



FOR MORE INFORMATION WWW.BEACHWHEELCHAIR.COM

DeBug Mobility Products
78 E Nine Mile Road
Pensacola, FL, 32534

850-478-5765
office@debugmobilityproducts.com
www.beachwheelchair.com

DEBUG XL FRAME

ALL-TERRAIN WHEELCHAIR

The Debug XL Beach Wheelchair features a reinforced frame engineered to support higher weight capacities, enhancing both stability and user confidence. Its unique five-wheel design evenly disperses weight, making it easier to push across soft sand. The seat is constructed from Phifertex material, triple-stitched with Tenera thread for exceptional durability in harsh outdoor conditions. Together, these features make the Debug XL a reliable choice for beach mobility, combining structural integrity with ease of use.

- LONG LASTING #316L STAINLESS STEEL FRAME
- REMOVABLE TELESCOPIC FOOTREST
- FITS THROUGH A STANDARD ADA DOORWAY
- POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS GRASS, GRAVEL, SAND, AND SNOW
- TRIPLE STITCHED SLING WITH TENERA 130 THREAD
- PHIFERTEX PLUS SLING MATERIAL IS MACHINE WASHABLE AND RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE

SPECIFICATIONS

Seat Sizes	ENGLISH (IN)	METRIC (CM)
Widths	20, 22	51, 56
Depth	17	43
Seat-to-Floor height	17, 21	43, 53
Backrest Height	18	46
Telescopic Footrest <i>Knee-to-Heel</i>	15-22	38-56
Weight Capacity	350-500 lbs	159-227 kg
Chair Width	33.5 inches	85 cm
Chair Weight	70-75 lbs	32-34 kg