

DEBUG ELEVATING LEG REST ALL-TERRAIN WHEELCHAIR



FOR MORE INFORMATION WWW.BEACHWHEELCHAIR.COM

DeBug Mobility Products
78 E Nine Mile Road
Pensacola, FL, 32534

850-478-5765
office@debugmobilityproducts.com
www.beachwheelchair.com

DEBUG ELEVATING LEG REST

ALL-TERRAIN WHEELCHAIR

The DeBug Beach Wheelchair delivers smooth, effortless mobility on sand, grass, and uneven terrain with its signature Wheeleez™ balloon wheels. The elevating leg rest provides crucial support for users with circulatory challenges, offering improved comfort, reduced pressure, and relief from leg spasms. Built for stability, ease of use, and true beach accessibility, it's the chair that lets everyone enjoy the shoreline in comfort.

- LONG LASTING #316L STAINLESS STEEL FRAME
- REMOVABLE TELESCOPIC FOOTREST
- ARTICULATING LEGREST
- FITS THROUGH A STANDARD ADA DOORWAY
- POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS GRASS, GRAVEL, SAND, AND SNOW
- DUAL STITCHED SLING WITH TENERA 130 THREAD
- PHIFERTEX PLUS SLING MATERIAL IS MACHINE WASHABLE AND RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE

SPECIFICATIONS

Seat Sizes	ENGLISH (IN)	METRIC (CM)
Widths	16, 18, 22	41, 46, 56
Depth	17	43
Seat-to-Floor height	17, 21	43, 53
Backrest Height	18	46
Telescopic Footrest <i>Knee-to-Heel</i>	15-22	38-56
Weight Capacity	Up to 300 lbs	Up to 136 kg
Chair Width	33.5 inches	85 cm
Chair Weight	60-65 lbs	27-29 kg