

DEBUG BEACH ROLLATOR



FOR MORE INFORMATION WWW.BEACHWHEELCHAIR.COM

DeBug Mobility Products
78 E Nine Mile Road
Pensacola, FL, 32534

850-478-5765
office@debugmobilityproducts.com
www.beachwheelchair.com

DEBUG BEACH ROLLATOR

The Debug Aluminum Beach Walker provides stability for individuals while walking on sandy beaches, rough terrain, and that normal walkers cannot accommodate. Light weight, portable, and ready to go when you are.

- LIGHTWEIGHT ALUMINUM FRAME FOLDS FOR TRANSPORTATION AND STORAGE
- HAND BRAKE AND LOCKING MECHANISM
- STAINLESS STEEL CASTERS AND AXLES RESIST SALT AND SAND
- PADDED SEAT AND BACKREST / UNDER SEAT STORAGE POUCH
- PNEUMATIC WHEELS ROLL ACROSS SAND, SNOW, AND GRAVEL
- LOW MAINTENANCE DELRIN BUSHINGS
- ADJUSTABLE AND REMOVABLE PUSH HANDLES

SPECIFICATIONS

	ENGLISH (IN)	METRIC (CM)
Handle Height	36 - 43	91 -109
Overall Width	39	99
Seat height	24 - 27	61 - 69
Weight	40lbs	18kg
Weight Capacity	Up to 350 lbs	Up to 136 kg

DeBug Beach Rollator Walker is intended for outside use only. It may not easily fit through doorways. While it is transportable and folds, it is also heavier than standard rollators. Some users may require assistance with setup and transportation .