

DEBUG FIXED LEG REST ALL-TERRAIN WHEELCHAIR



FOR MORE INFORMATION WWW.BEACHWHEELCHAIR.COM

DeBug Mobility Products
78 E Nine Mile Road
Pensacola, FL, 32534

850-478-5765
office@debugmobilityproducts.com
www.beachwheelchair.com

DEBUG FIXED LEG REST ALL-TERRAIN WHEELCHAIR

The DeBug Fixed Leg Rest Beach Wheelchair is designed to provide dependable mobility across terrain where conventional wheelchairs cannot travel. Whether exploring sandy beaches, snow-covered sidewalks, grassy parks, or forest trails, the DeBug delivers a smooth and stable ride. Constructed with durable Phifertex sling seating, dual-stitched with Tenera thread, it's built to withstand years of outdoor use.

- LONG LASTING #316L STAINLESS STEEL FRAME
- REMOVABLE TELESCOPIC FOOTREST
- FITS THROUGH A STANDARD ADA DOORWAY
- POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS GRASS, GRAVEL, SAND, OR SNOW
- DUAL STITCHED SLING WITH TENERA 130 THREAD
- PHIFERTEX PLUS SLING MATERIAL IS MACHINE WASHABLE AND RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE
- CUSTOM BUILT TO FIT CHAIR USER'S MEASUREMENTS

SPECIFICATIONS

Seat Sizes	ENGLISH (IN)	METRIC (CM)
Widths	16, 18, 22	41, 46, 56
Depth	17	43
Seat-to-Floor height	17, 21	43, 53
Backrest Height	18	46
Telescopic Footrest <i>Knee-to-Heel</i>	15-22	38-56
Weight Capacity	Up to 300 lbs	Up to 136 kg
Chair Width	33.5 inches	85 cm
Chair Weight	60-65lbs	27-29 kg