

DEBUG JUNIOR FIXED LEG REST ALL-TERRAIN WHEELCHAIR



FOR MORE INFORMATION WWW.BEACHWHEELCHAIR.COM

DeBug Mobility Products
78 E Nine Mile Road
Pensacola, FL, 32534

850-478-5765
office@debugmobilityproducts.com
www.beachwheelchair.com

DEBUG JUNIOR (FIXED LEG REST) ALL-TERRAIN WHEELCHAIR

The DeBug Junior Beach Wheelchair is designed specifically for children and adolescents who need reliable mobility across sand and other soft terrain. With a smaller frame and seating dimensions, the Junior model provides the same durability and performance as our full-size beach wheelchairs while offering a better fit for younger users. Every chair is custom built to order and can be tailored to allow room for growth as your child develops.

- LONG LASTING #316L STAINLESS STEEL FRAME
- REMOVABLE TELESCOPIC FOOTREST
- FITS THROUGH A STANDARD ADA DOORWAY
- POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS GRASS, GRAVEL, SAND, AND SNOW
- DUAL STITCHED SLING WITH TENERA 130 THREAD
- PHIFERTEX PLUS SLING MATERIAL IS MACHINE WASHABLE AND RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE

SPECIFICATIONS

Seat Sizes	ENGLISH (IN)	METRIC (CM)
Width	16	41
Depth	17	43
Seat-to-Floor height	17, 21	43, 53
Backrest Height	18	46
Telescopic Footrest <i>Knee-to-Heel</i>	15-22	38-56
Weight Capacity	Up to 250 lbs	Up to 113 kg
Chair Width	33.5 inches	85 cm
Chair Weight	60-65lbs	27-29 kg